

Professional results for a fraction of the price!



*"About time! Finally a product that does what it says it can do! Yeah!" - Michelle*

DERMANESSE™

2261 Brookhollow Plaza Drive, Suite 200  
Arlington, TX 76006

Customer Service 866-849-5998  
Order Phone 800-379-3800  
[www.FabulousSkin.com](http://www.FabulousSkin.com)

\*All quotes are from actual testimonials on file. DermaSponge is Patent Pending.

DERMA  SPONGE  
A Younger You™



**The First Professional  
Microdermabrasion  
System for Home Use!  
Results Guaranteed!**

## Congratulations!



*“Since using DermaSponge my skin looks terrific!” - Susan*

## You are about to give yourself a professional Microdermabrasion Treatment!

The six treatment DermaSponge home microdermabrasion system includes everything you need to achieve optimum results for much less than just one treatment at a doctor’s office, spa or salon. DermaSponge utilizes the same aluminum oxide crystals used in professional microdermabrasion systems—for professional results in the convenience of your own home!

DermaSponge is safe for all skin types and works great for either men or women. It’s fast and easy, too. In **just 3 minutes every two weeks**, DermaSponge will help brighten your skin, minimize fine lines and wrinkles, reduce large pores, help acne scars disappear and even-out skin tone.

In many cases, you’ll see results with your first treatment. However, since everyone’s skin is different, individual results may vary.



## **Important!** Please read before you begin.

**M**icrodermabrasion is a highly effective method of improving your skin's appearance and health. By gently removing the surface layer of skin, DermaSponge stimulates growth of new, younger looking, more radiant skin. Until DermaSponge, microdermabrasion was only available in doctors' offices or clinics. Now you can achieve the same results at home!

***To get the best results, you must read all the instructions thoroughly and follow them carefully! Failure to follow the instructions may result in ineffective treatments or excessive skin irritation! Please read the instructions before each treatment.***

Since everyone's skin varies, we recommend starting with a 1 minute application using light to moderate pressure. For unusually sensitive skin, perform the Skin Sensitivity Test on the next page first. As your skin adjusts, you may extend your time a little each treatment until you reach 3 minutes.

If you have little or no irritation, redness, or peeling after a test or treatment, your skin may be less sensitive or may be adjusting to the treatments. If that happens, increase treatment times gradually up to 5 minutes over the entire face. Always start with lighter pressure and adjust for your skin type as needed. Thicker, less sensitive skin may be able to handle slightly more pressure.

***"I was so surprised how great it worked, I cancelled an expensive microdermabrasion appointment at the spa. Thank you for an affordable alternative!" - Lori***

### ***Skin Sensitivity Test***

*Use DermaSponge on the soft skin of your inner forearm for a minute. Wait awhile to see the results. If you have irritation and redness similar to a minor sunburn, this is normal after true microdermabrasion. If you experience more discomfort, then you may have sensitive skin and should limit your first couple of treatments to 1 minute. For extremely sensitive skin, apply HydraLipid Serum before sponge application to reduce the irritation.*

***Important:*** Please do not begin DermaSponge treatments if you have an open wound, infection, sores, inflamed acne, or other skin irritations. For optimal results allow your skin to heal completely before using DermaSponge.

**Note:** there are no known allergic reactions to aluminum oxide crystals. No harsh chemicals are ever used. Do not use abrasive products, such as Retin-A®, Renova® or Alpha Hydroxy Acids on the skin for several days after using DermaSponge.



## Skin Preparation

*Apply DermaSponge Prep Wash Cleanser.* DermaSponge Prep Wash is designed to cleanse your skin before microdermabrasion with gentle, soothing ingredients like Eucalyptus Oil, Grapefruit Seed Extract, Sandalwood Extract, Lemon Oil and Apricots. Avoid eye area. Does not lather. Thoroughly dry your skin.

*For extremely sensitive skin only:* Apply HydraLipid Serum before sponge application to reduce irritation.



*“I want to use this product forever. Forget all those chemicals that promise youthfulness—this really works!” - Donna*

## Sponge Application

Use dry sponge only on dry skin. Use one sponge per treatment. *Do not reuse sponge. Discard after use for proper hygiene.*

Following the examples on the next page, gently buff skin with the abrasive side of the sponge using light to moderate pressure to test the sensitivity of your skin.

If your skin is unusually sensitive and you experience excessive redness, burning or irritation, limit your treatment to 1 minute. As you do more treatments, your skin may adjust, allowing you to gradually increase application times.

If there is little or no irritation, you may increase application time. Do not go beyond 5 minutes for your whole face as it is not necessary and may irritate the skin excessively.

*Never use heavy pressure or spend too much time on just one area of the skin.*

**Around Eyes & Lip Line** - buff very lightly with extreme caution as the skin around the eye area is very thin. Use extremely light, gentle circular motions, buffing away from the center of the face outwards.

Remember, DermaSponge can also be used on the neck, chest and the hands as well as other areas of the body.



# How to do a Derma Sponge Treatment



## 1. Forehead

Using a circular motion, buff your skin in small circles moving continuously from one side to the other across your forehead.

## 3. Under the Nose

Buff horizontally back and forth between your nose and upper lip.



## 2. Cheeks

Buff in a broad circular motion from the front towards the back on both cheeks.



## 4. Chin

Buff horizontally back and forth across your chin.

## After Your DermaSponge Treatment



1. Apply DermaRelief Gel generously all over your face. Derma Relief Gel contains Aloe Vera and lidocaine to soothe skin after your treatment. Apply as often as needed to relieve irritation. Avoid eye area.

2. Apply HydraLipid Serum to hydrate your skin. This exclusive product restores and hydrates your skin after microdermabrasion and with daily use. You can also use HydraLipid Serum at night on cleansed skin. Use alone or under a rich moisturizer for extra protection and hydration. Avoid eye area.



3. Apply Protective Day Creme SPF 15 every day to protect your fresh, new skin from damaging sun exposure. It is also a wonderful base for under make-up. This rich moisturizer contains Aloe Vera Gel, Vitamins E and A, and Green Tea Extracts to protect your skin from free radical and sun damage. Avoid eye area.



Note: About 2-3 days after treatment, your skin may start to peel or flake like after a mild sunburn. This is a normal healing process for some people after real microdermabrasion and will dissipate within a day or so. If this happens, just gently rub the dead skin away with a towel using warm water.



*"You changed my life and saved me tons of money! I have been getting compliments like crazy! Thank you for giving me my confidence back!" - Carmel*

*Unconditional Guarantee - Treat yourself to two treatments on us and see the results for yourself. If you are not completely satisfied for any reason, just return the complete kit less any used sponges within 30 days for a full refund (less shipping). No questions asked!*

## Frequently Asked Questions

### *What is microdermabrasion?*

Originally developed in Italy, one of the beauty capitals of the world, microdermabrasion has been used successfully for years by millions of women for the treatment of fine lines, wrinkles, acne scars, superficial pigmentation and sun-damaged skin. It is considered safe on all skin types. Though professional treatments can cost thousands of dollars, DermaSponge is the only home product that gives you results comparable to professional treatments for a fraction of the cost.

### *Will DermaSponge cause redness on my skin?*

Immediately after treatment, your skin will take on a fresh new "glow" and there will be some temporary redness present, similar to a mild sunburn. This will disappear in a few hours (or days for more sensitive skin). Apply DermaRelief Gel as often as necessary to relieve irritation. Makeup can be applied immediately after procedure to minimize the appearance of redness.

### *Can I still use Retin-A®, Renova® and Alpha Hydroxy Acids?*

Yes. But DO NOT USE such abrasive products or chemical exfoliants on your skin for several days after using DermaSponge.

### *Why does my skin flake or peel after using DermaSponge?*

This is good! Your new, younger skin is healing and rejuvenating. Like after a mild sunburn, the peeling and flaking should last about 2-3 days (a little more for some people with sensitive skin). Just gently rub the dead skin away with a soft washcloth or loofa and warm water.

### *Will DermaSponge increase my skin's sensitivity to the sun?*

Not directly, since there are no chemicals in DermaSponge. You should use a sunscreen for several days after any true microdermabrasion treatment as your new, younger-looking skin will be slightly more sensitive than normal.

### *How much pressure should I use?*

Never use heavy pressure. Use gentle, circular motions over the entire face. If you have sensitive skin, test DermaSponge on the underside of your forearm for a minute first. If you have excessive redness and irritation, use lighter pressure. As your skin adjusts you may be able to gradually increase the pressure and time, up to 3 minutes, or 5 minutes for tougher skin.

### *My skin is burning or tingling - is this OK?*

Yes. As with any true microdermabrasion treatment, you may experience redness, tingling or burning sensations similar to a mild sunburn. This means your microdermabrasion treatment is working! This will go away in a few hours. Apply DermaRelief Gel as often as necessary. With more sensitive skin, it can take a little longer. Remember, if your skin is extremely sensitive, you can use the HydraLipid Serum *before* your DermaSponge treatment to reduce irritation.

*“Just one treatment gave me the same results as a professional microdermabrasion treatment! Thanks!” - Allison*

**Everything you need to keep your new skin looking young!**  
**Join our "Younger You" Club now and get FREE Shipping!**

**Deluxe DermaSponge Kit** - includes 6 DermaSponges embedded with the same aluminum oxide crystals used by professionals, plus DermaRelief Gel, our mild Prep Wash Cleanser, rich Hydralipid Serum moisturizer, and Protective Day Crème SPF 15 sunscreen, all optimized for use with our microdermabrasion treatments.

**DermaRelief Gel** (1.7 oz./50 ml) With Aloe Vera Gel and Lidocaine for the ultimate in soothing the skin.

**Protective Day Crème SPF 15** (1.7 oz./50 ml) This moisturizer contains Aloe Vera Gel, Vitamins E and A, and Green Tea Extracts.

**PrepWash Cleanser** (1.7 oz./50 ml) A soothing, gentle cleanser with Eucalyptus Oil, Grapefruit Seed Extract, Sandalwood Extract, Lemon Oil, and Apricots.

**HydraLipid Serum** (1/2 oz.) To restore and hydrate the skin after treatment and every day usage.

**Join our "Younger You" Club & get FREE Shipping** (USA & Canada; or just \$6.95 International). We'll send you an AutoShip every 90 days so you can always keep your new skin looking young. Cancel at any time online or call our toll-free phone number.

Customer Service 866-849-5998  
 Order Phone 800-379-3800

**Order Form**

Item	Price	Quantity	Subtotal
<b>Deluxe Microdermabrasion Kit</b> <i>Contains 6 sponges, plus 1 each of the following:</i>	<b>\$79.95</b>		
<b>DermaRelief Gel</b>	<b>\$24.00</b>		
<b>Prep Wash Cleanser</b>	<b>\$27.00</b>		
<b>Hydralipid Serum</b>	<b>\$30.00</b>		
<b>Protective Day Crème SPF 15</b>	<b>\$33.00</b>		
<b>Subtotal</b>			

**YES.** I want to join the "Younger You" Club and get FREE Shipping USA & Canada; or just \$6.95 International. I'll get one kit every 90 days, starting 3 months after this order. I may cancel at any time.

**FL Residents** Add 6% Sales Tax \_\_\_\_\_  
**Shipping:** \$5.95 USA & Canada \_\_\_\_\_  
 \$12.95 International \_\_\_\_\_  
**TOTAL** \_\_\_\_\_

First Name: \_\_\_\_\_ Last Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State/Province: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Email Address: \_\_\_\_\_ Daytime Phone: \_\_\_\_\_ Nighttime Phone: \_\_\_\_\_

Payment Type:  Visa  Mastercard  Amex  Discover  Check/Money Order

Name on Card: \_\_\_\_\_

Billing address (if different than above):  
 \_\_\_\_\_

Credit Card #: \_\_\_\_\_ Expiration Date: \_\_\_\_\_

Signature \_\_\_\_\_